

# RQI Cleaning, Decontamination And Hand Hygiene

Materials such as adult and infant bag-valve masks (BVM), manikin cleaning wipes (PDI Super Sani-Cloth), hand sanitizer, and other consumables should be placed in the storage bin located under the adult manikin platform. Any non-functioning components should be reported as soon as observed for repair or replacement as required.

**Writing utensils of any kind are strictly prohibited around the RQI cart.**

1. Before using the RQI platform hand sanitizer must be used. \*
2. After completion gloves are to be worn during cleaning of the RQI Station. \*
3. Manikins, computer keyboard and any hard surface touched will need to be cleaned with PDI Super Sani-Cloths (Purple Top) cleaning wipes which are stored in the storage bin. Wipe any hard surface that was touched, make sure to clean the chest and the mouth of both adult and infant manikin. Allow 2 minutes for the surfaces to dry before use again. Gloves must be worn during this cleaning. \* (Supply Chain Update attached for Clorox Wipes use.) \*\*\*
4. All materials removed from the storage bin should now be replaced and the bin put back underneath the adult manikin.
5. After cleaning has been completed all wipes need to be put in the general trash can with your gloves which were removed per IU Health protocol. These do not need to go into Biohazard (Red) Bags, just general trash.
6. Wash your hands after completion of all of these steps per IU Health Protocol. \*\*

\*- IU Health Policy Stat ID 6161723

\*\* - IU Health Policy Stat ID 5624388

\*\*\* Page 2



# RQI Supply Chain Updates

## Disinfectant Wipes

<b>Situation</b>	Purple and Grey top wipes are in low supply.
<b>Background</b>	RQI Manikin cleaning and decontamination practices are an essential part of infection prevention practice. The current recommendation includes the use of Purple or Grey top wipes in this practice.
<b>Assessment</b>	An alternate resource is needed to replace Purple and Grey Top wipes.
<b>Recommendations</b>	Please, use Clorox wipes as pictured below. Do not use Bleach containing wipes on any electronics.

### Hot Topics 2: Supply Chain Updates



**Sporicidal Disinfectant**  
Use for  
**Enteric Isolation/C. diff**



Contains Bleach

3 Minute Contact Time

**General Purpose Disinfectant –**  
Substitute for Purple Top and Grey To  
wipes



Bleach Free

4 Minute Contact Time

3:50

