

Fuji Synapse PACS Quick Guide for Providers

Use your IU Health User Name and Password to access Fuji Synapse PACS

Need your password reset?
Call 317-962-2828 or your local helpdesk number


4. Select the patient

Enter the patient's MRN then click 'Enter'

Patient Name	Patient ID
	123456789

1. Open Synapse

From your desktop, double click on the Synapse icon to view the available databases.



Double click
IU Health

5. Viewing Images

Open the study: Double click on the study to view the images.


Scrolling through images one at a time:

Left click on the image and scroll the scroll wheel

Scrolling through images quickly:

Press down/hold the scroll wheel and move the mouse back and forth.

2. Log in



6. Power Jacket

To open the Power Jacket: Click the Power Jacket icon

To view study info: Click the Study Info icon

To view the report: Click the Show Report Icon


To view or add a note: Click the Show Notes icon

To view the series: Click the Series Thumbnails icon

7. Viewing Comparisons: From Power Jacket, double click on another exam and it will display as a comparison OR select the comparison exam and then click Compare

Compare

3. Select Folder




'All Recent Studies with Images' displays studies with images from past 7 days.

'All Studies' displays studies more than 7 days old

8. Image Tools

Image tools include Window Levels, Zoom, Ruler, Text and more

To access Image Tools: Rt click on an image



Accessing Images from Cerner PowerChart:

The 'Image Link' icon is available in PowerChart for Radiology Studies with images.

- It can be found on the toolbar once the Radiology exam is opened from Results Review, Clinical Notes or the ED Tracking Shell.
- Double clicking the 'Image Link' icon will load the images associated to the selected exam.

